

PRODUCT CATALOGUE

SepaSperm®

- For sperm preparation by the density gradient centrifugation.
- Free from animal derivatives.
- SepaSperm® Wash Solution is HEPES buffered.



	REF	Code	Contents
SepaSperm® Solution 100%	92143	SE100G-100	100mL
SepaSperm® Solution 80%	92073	SE80G-100	100mL
	92079	SE80G-50	50mL
SepaSperm® Solution 60%	92074	SE60G-100	100mL
	92080	SE60G-50	50mL
SepaSperm® Solution 40%	92075	SE40G-100	100mL
	92081	SE40G-50	50mL
SepaSperm® Wash Solution	92147	SEWSG-100	100mL

COMPONENTS



Calcium chloride / Gentamicin sulfate / Glucose / HEPES /
Magnesium sulfate / Potassium chloride / Potassium phosphate /
Recombinant human albumin / Silica particles / Sodium bicarbonate /
Sodium chloride / Sodium lactate / Sodium pyruvate

QUALITY CONTROL

pH 7.2-7.6 / Osmolarity 270-300 mOsm/L / Endotoxin <0.25EU/mL /
Mouse Embryo Assay ≥80% / Sperm Survival (24h) ≥80% /
Sperm Penetration ≥3 / Sterility Test
Storage: 2-8°C
Shelf life: 18 months

RESULT

Results of SepaSperm®



	N	Sperm Concentration (×10 ⁶ /mL)	Sperm motility rate (%)	Sperm motility recovery rate (%)
Semen	5	128.7±40.0	58.7±19.7	—
Prepared with SepaSperm® Solution	5	126.5±50.4	91.5±6.5	26.6±7.8
Prepared with Product A	5	114.7±55.1	90.3±7.9	22.4±5.5

Numbers in mean±SD

Reference: Kusuha Women's Clinic

Specification may change without pre-notice for purpose of product improvement.

Kitazato Corporation


HEADQUARTERS : 100-10 Yanagishima, Fuji, Shizuoka 416-0932 JAPAN

TOKYO : 1-1-8 Shibadaimon, Minato-ku, Tokyo 105-0012 JAPAN

Mail contact@kitazato.co.jp

An example of single-layer gradient centrifugation

- 1 Bring all solutions to room temperature (20-25°C) before use. Prepare two 15mL centrifuge tubes. One of them with 3-5mL of SepaSperm® 80%. **A** The other with 5mL of SepaSperm® Wash Solution. **B**
- 2 Check the semen volume with a 5mL pipette or similar and evaluate the semen sample with a sperm count chamber.
- 3 Pipette semen on top of the prepared SepaSperm® 80%. Stir 2cm vertically around the layers between the semen and SepaSperm® 80%.
- 4 Centrifuge for 15 minutes at approximately 600 x g. (Separation of sperm and seminal fluid)
Depending on the state of collected sperm, adjust xg or centrifugation time accordingly.
- 5 After centrifugation, remove the supernatant with a new pipette.
- 6 Add 5mL of SepaSperm® Wash Solution (SEWS) to the sperm pellet and resuspend. Centrifuge at 400 x g for 5 minutes for washing.
Depending on the state of collected sperm, adjust xg or centrifugation time accordingly.
- 7 Remove the supernatant leaving the sperm pellet only. Dilute the sperm pellet with SepaSperm® Wash Solution according to the following treatment.



An example of two-layer density gradient centrifugation

- 1 Bring all the solutions to room temperature (20-25°C) before use. Prepare two 15mL centrifuge tubes. One of them with layers in the next order from bottom to up ;
 - 1 SepaSperm® 80% 1.5-2.0mL
 - 2 SepaSperm® 40% 1.5-2.0mL
 - 3 Liquefied semen
The other with 5mL of SepaSperm® Wash Solution.
- 2 Check the semen volume with a 5mL pipette or similar and evaluate the semen sample with a sperm count chamber.
- 3 Centrifuge the tube at 600 x g for 15 minutes. (Separation of sperm and seminal fluid)
Depending on the state of collected sperm, adjust xg or centrifugation time accordingly.
- 4 Aspirate the sperm pellet.
- 5 Transfer the sperm pellet to a new centrifuge tube prepared with 5mL of SepaSperm® Wash Solution, and mix well.
- 6 Centrifuge at 400 x g for 5 minutes for washing.
Depending on the state of collected sperm, adjust xg or centrifugation time accordingly.
- 7 Remove the supernatant leaving the sperm pellet only. Dilute the sperm pellet with SepaSperm® Wash Solution according to the following treatment.

