


PRODUCT CATALOGUE

# Oocyte Cryo Plate

- Premarkings on the plate guide you to prepare the different solutions in appropriate volumes.
- High repeatability in the equilibration conditions leads to stable results.
- Ditches on the lid hold Cryotop stable for better handling.



**Collaborative Development** : Kato Ladies Clinic / Keio University, Department of. Obstetrics & Gynecology / National Center for Child Health and Development / Obstetrics and Gynecology, St Marianna University School of Medicine / Obstetrics and Gynecology, Graduate School of Medicine and Faculty of Medicine, The University of Tokyo


REF	Code	Contents
83061	Oocyte Cryo Plate	10pcs/pk

## QUALITY CONTROL

Sterility Test / Endotoxin  $\leq 0.5\text{EU/unit(EU/mL)}$  / Mouse Embryo Assay  $\geq 80\%$

## RESEARCH DATA

Comparison of osmolality change between the conventional and the new protocol with the Oocyte Cryo Plate.



- Oocyte Cryo Plate allows solutions to mix moderately. It allows oocytes to equilibrate more gradually compared to the conventional protocol.
- Gradual osmolality change improve oocyte development after warming.

Specification may change without pre-notice for purpose of product improvement.

## Kitazato Corporation

HEADQUARTERS : 100-10 Yanagishima, Fuji, Shizuoka 416-0932 JAPAN

TOKYO : 1-1-8 Shibadaimon, Minato-ku, Tokyo 105-0012 JAPAN

**Mail** [contact@kitazato.co.jp](mailto:contact@kitazato.co.jp)

# PROTOCOL

The video protocol is available  
on our official YouTube channel.




Fig.1




Fig.2




Fig.3





Fig.4




Fig.5




Fig.6

Perform the following pre-freezing  
procedures using our oocyte/embryo  
vitrification solution.

- 01** On Oocytes Cryo Plate, dispense 20  $\mu$ L drop of BS in area ①, 300  $\mu$ L of ES in area ②, 300  $\mu$ L of VS in area ③ and ④ respectively. Keep the lid on until use (Fig.1).
- 02** Transfer the oocytes from the culture medium to the BS drop.
- 03** Dispense 20  $\mu$ L of ES drop on the right side of the BS drop, and merge the ES to the BS drop using the chip of the pipette, then leave it still for 3 min (Fig. 2).
- 04** Dispense 20  $\mu$ L drop of ES below the mixed drop. Merge the new ES drop to it with the chip, and leave it still for 3 min (Fig. 3).
- 05** Transfer the oocytes to the surface of the dispensed ES in area ② and leave it still for 9 min (Fig. 4).
- 06** After completion of ES equilibration, aspirate the oocytes/embryos with minimal volume of ES at the tip of the pipette . Transfer the aspirated oocytes /embryos to the surface of VS1 and set the timer.  
\*Referential operation time in VS1 and VS2 together is 60 - 90 seconds to ensure dehydration. It is recommended to keep oocytes/embryos immersed in VS for at least 60 - 90 seconds.
- 07** Expel the remaining ES in the pipette to the outside the well, and wash the pipette by aspirating and expelling sufficient volume of VS.
- 08** Aspirate the oocytes/embryos and perform steps the below ① and ② in three positions in the VS1.
  - ① Expel the oocytes/embryos from the pipette to the VS.
  - ② Stir around the oocytes/embryos 5 times gently for 5 seconds.
 The oocytes/embryos are equilibrated and dehydrated by performing steps ① and ② in three different positions (Fig. 5).

- 09** Wash the pipette with VS2 as explained in STEP2. Perform ① and ② of step 08 twice at different positions in the VS2. Confirm the completion of VS equilibration by the below two criteria before loading oocytes/embryos on Cryotop.
  - The oocytes/embryos are shrunk.
  - The oocytes/embryos does not rise to the surface. (They stay in focus under microscope.)

If the above two points are not confirmed, perform ① and ② again in no rush (Fig. 6).