

PRODUCT CATALOGUE

Flushing Buffer

- For flushing follicles.
- To be used in ambient atmosphere for washing oocytes in the retrieval.
- Flushing Buffer does not contain protein to avoid foaming.



Order Number	Code	Contents
93451	ASP-100	100mL
93450	ASP-500	500mL

COMPONENTS

Calcium chloride / Essential amino acids / Gentamicin /
Glucose / Heparin / HEPES / Magnesium sulfate /
Non-essential amino acids / Potassium chloride /
Potassium phosphate / Sodium bicarbonate /
Sodium chloride / Sodium lactate / Sodium pyruvate /
Taurine

QUALITY CONTROL

pH 7.2-7.6 / Osmolality 270-295mOsm/kg /
Endotoxin <0.25EU/mL / Sterile Filtration / Sterility Test /
Mouse Embryo Assay ≥80%
Storage: 2-8°C
Shelf Life: 180 days

Specification may change without pre-notice for purpose of product improvement.

Kitazato Corporation

HEADQUARTERS : 100-10 Yanagishima, Fuji, Shizuoka 416-0932 JAPAN

TOKYO : 1-1-8 Shibadaimon, Minato-ku, Tokyo 105-0012 JAPAN

Mail contact@kitazato.co.jp