

## PRODUCT CATALOGUE

## Hypure® Oil

- Light and Heavy, 2 types are available.
- High-purity paraffin oil without RNA detection.
- Prepared oil which is washed twice with ultrapure water.
- Light shielding glass is used to prevent deterioration of quality by light. Air tightness is kept in the bottle container.



## Hypure® Oil Heavy

Order Number	Code	Contents
93620	MOHV-500	500mL
93621	MOHV-100	100mL

## Hypure® Oil Light

Order Number	Code	Contents
93520	MOLH-500	500mL
93521	MOLH-100	100mL

## COMPONENTS

High-purity Paraffin Oil

Specific Gravity

Heavy : 0.862g/mL (25°C) Light : 0.838g/mL (25°C)

## QUALITY CONTROL

Endotoxin <0.25EU/mL / Sterility test /

Mouse Embryo Assay ≥80% / POV ≤0.1meq/kg /

Readily Carbonizable Substances/

Storage: 8-30°C

Shelf life: 18 months

Specification may change without pre-notice for purpose of product improvement.

## Kitazato Corporation

HEADQUARTERS : 100-10 Yanagishima, Fuji, Shizuoka 416-0932 JAPAN

TOKYO : 1-1-8 Shibadaimon, Minato-ku, Tokyo 105-0012 JAPAN

[Mail contact@kitazato.co.jp](mailto:contact@kitazato.co.jp)