

PRODUCT CATALOGUE

PBS(-)
Phosphate buffered saline pH7.4

- Without calcium and magnesium.
- Suitable for washing cells, especially for washing oocytes.
- Suitable to use under ambient atmosphere.



REF	Code	Contents
93352	PBSF-500	500mL
93353	PBSF-100	100mL
93354	PBSFH-500	500mL
93355	PBSFH-100	100mL

*PBSFH is with heparin.

COMPONENTS

Di-sodium hydrogen phosphate anhydrous /
Potassium chloride / Potassium phosphate /
Sodium chloride
PBSFH: Heparin

QUALITY CONTROL

pH 7.2-7.6 / Osmolarity 270-295mOsm/L /
Endotoxin <0.25EU/mL /
Mouse Embryo Assay ≥80% / Sterility Test
Storage: 15-30°C
Shelf Life: 2 years

Specification may change without pre-notice for purpose of product improvement.

Kitazato Corporation

HEADQUARTERS : 100-10 Yanagishima, Fuji, Shizuoka 416-0932 JAPAN
TOKYO : 1-1-8 Shibadaimon, Minato-ku, Tokyo 105-0012 JAPAN

Mail contact@kitazato.co.jp